

Product Datasheet

Mouse PYY ELISA Kit (orb777223)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	Q9EPS2
Sensitivity	4.52 pg/mL
Expiration Date	Please enquire.

