

## Product Datasheet

### Mouse ACE2 ELISA Kit (orb777589)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Mouse
<b>Range</b>	12.5-800 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Uniprot ID</b>	<b>Q8R0I0</b>
<b>Sensitivity</b>	5.5 pg/mL
<b>Expiration Date</b>	Please enquire.

