

Product Datasheet

Human IL17D ELISA Kit (orb777767)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	7.82-500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q8TAD2
Sensitivity	2.7 pg/mL
Expiration Date	Please enquire.

