

## Product Datasheet

### Human APLN ELISA Kit (orb777933)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Human
<b>Range</b>	25-1600 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma and other biological fluids
<b>Uniprot ID</b>	<b>Q9ULZ1</b>
<b>Sensitivity</b>	8.25 pg/mL
<b>Expiration Date</b>	Please enquire.

