

## Product Datasheet

### Porcine HB ELISA Kit (orb778122)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Porcine
<b>Range</b>	4.69-300 µg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, erythrocyte lysates and other biological fluids
<b>Sensitivity</b>	1.74 µg/mL
<b>Expiration Date</b>	Please enquire.

