

Product Datasheet

Mouse CD163 ELISA Kit (orb778347)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Mouse |
| Range | 31.25-2000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma and other biological fluids |
| Uniprot ID | Q2VLH6 |
| Sensitivity | 13.3 pg/mL |
| Expiration Date | Please enquire. |

