

## Product Datasheet

### Gallus GHRL ELISA Kit (orb778366)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Gallus
<b>Range</b>	156.25-10000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	46.4 pg/mL
<b>Expiration Date</b>	Please enquire.

