

Product Datasheet

Mouse TLR6 ELISA Kit (orb778558)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Mouse |
| Range | 1.57-100 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Tissue homogenates, cell lysates, cell culture supernates and other biological fluids. |
| Uniprot ID | Q9EPW9 |
| Sensitivity | 0.54 ng/mL |
| Expiration Date | Please enquire. |

