

## Product Datasheet

### Human ABCC8 ELISA Kit (orb778592)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 0.32-20 ng/mL   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Sample Types</b>        | tissue homogenates, cell lysates and other biological fluids  |
| <b>Uniprot ID</b>          | <a href="#">Q09428</a>  |
| <b>Sensitivity</b>         | 0.52 ng/mL  |
| <b>Expiration Date</b>     | Please enquire.   |

