

Product Datasheet

Human PLA2G5 ELISA Kit (orb778640)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 78.13-5000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids |
| Uniprot ID | P39877 |
| Sensitivity | 31 pg/mL |
| Expiration Date | Please enquire. |

