

Product Datasheet

Human FLRT3 ELISA Kit (orb778797)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 1.57-100 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Tissue homogenates, cell lysates and other biological fluids. |
| Uniprot ID | Q9NZU0 |
| Sensitivity | 0.57 ng/mL |
| Expiration Date | Please enquire. |

