

Product Datasheet

Mouse FGF2 ELISA Kit (orb778804)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P15655
Sensitivity	4.41 pg/mL
Expiration Date	Please enquire.

