

Product Datasheet

Rat TTR ELISA Kit (orb778854)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P02767
Sensitivity	0.146 ng/mL
Expiration Date	Please enquire.

