

Product Datasheet

Human ADM2 ELISA Kit (orb778874)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Human
Range	312.5-20000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q7Z4H4
Sensitivity	91.5 pg/mL
Expiration Date	Please enquire.

