

Product Datasheet

Human TAC1 ELISA Kit (orb778892)

| | |
|----------------------------|---|
| Description | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| Reactivity | Human |
| Range | 156.25-10000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | P20366 |
| Sensitivity | 52.8 pg/mL |
| Expiration Date | Please enquire. |

