

Product Datasheet

Mouse TGFb2 ELISA Kit (orb779074)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids.
Uniprot ID	P27090
Sensitivity	8.93 pg/mL
Expiration Date	Please enquire.

