

## Product Datasheet

### Human ACTN3 ELISA Kit (orb779082)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	1.57-100 ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<b>Q08043</b>
<b>Sensitivity</b>	0.62 ng/mL
<b>Expiration Date</b>	Please enquire.

