

Product Datasheet

Rat FGF7 ELISA Kit (orb779092)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Rat
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q02195
Sensitivity	5.37 pg/mL
Expiration Date	Please enquire.

