

Product Datasheet

Rat FGF23 ELISA Kit (orb779108)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Rat
Range	3.13-200 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q8VI82
Sensitivity	0.92 pg/mL
Expiration Date	Please enquire.

