

Product Datasheet

Animal HSP70 ELISA Kit (orb779116)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Bovine
Range	3.13-200 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q3ZCH0
Sensitivity	1.27 ng/mL
Expiration Date	Please enquire.

