

Product Datasheet

Mouse MSTN ELISA Kit (orb779138)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	781.25-50000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	008689
Sensitivity	253.7 pg/mL
Expiration Date	Please enquire.

