

Product Datasheet

Mouse MEP1a ELISA Kit (orb779165)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	62.5-4000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P28825
Sensitivity	19.1 pg/mL
Expiration Date	Please enquire.

