

Product Datasheet

MET (phospho-Tyr1234/1235) antibody (FITC) (orb7792)

Description	MET (phospho-Tyr1234/1235) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MET around the phosphorylation site of Tyr1234/1235
Target	Phospho-Met (Tyr1234 + Tyr1235)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	153kDa
Uniprot ID	P08581
Expiration Date	12 months from date of receipt.