

Product Datasheet

Human HJV ELISA Kit (orb779266)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Human
Range	15.63-1000 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q6ZVN8
Sensitivity	4.93 ng/mL
Expiration Date	Please enquire.

