

Product Datasheet

Human HJV ELISA Kit (orb779288)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 31.25-2000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | Q6ZVN8 |
| Sensitivity | 12.7 pg/mL |
| Expiration Date | Please enquire. |

