

## Product Datasheet

### Animal TGFb1 ELISA Kit (orb779316)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Bovine
<b>Range</b>	15.63-1000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, platelet-poor plasma, tissue homogenates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<b>P18341</b>
<b>Sensitivity</b>	4.67 pg/mL
<b>Expiration Date</b>	Please enquire.

