

Product Datasheet

Rat MEP1a ELISA Kit (orb779367)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Rat
Range	62.5-4000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q64230
Sensitivity	15.2 pg/mL
Expiration Date	Please enquire.

