

## Product Datasheet

### Porcine NOS3 ELISA Kit (orb779460)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Porcine
<b>Range</b>	15.63-1000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<b>Q28969</b>
<b>Sensitivity</b>	6.1 pg/mL
<b>Expiration Date</b>	Please enquire.

