

## Product Datasheet

### Animal CRP ELISA Kit (orb779461)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Bovine
<b>Range</b>	6.25-400 ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	2.43 ng/mL
<b>Expiration Date</b>	Please enquire.

