

## Product Datasheet

### Porcine D2D ELISA Kit (orb779469)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Porcine
<b>Range</b>	78.13-5000 ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	plasma
<b>Sensitivity</b>	25.4 ng/mL
<b>Expiration Date</b>	Please enquire.

