

## Product Datasheet

### Canine EPO ELISA Kit (orb779508)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>          | Canine  |
| <b>Range</b>               | 31.25-2000 pg/mL  |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Sample Types</b>        | serum, plasma and other biological fluids   |
| <b>Uniprot ID</b>          | <b>P33707</b>   |
| <b>Sensitivity</b>         | 13.7 pg/mL  |
| <b>Expiration Date</b>     | Please enquire.   |

