

Product Datasheet

Porcine IL13 ELISA Kit (orb779533)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Porcine
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q95J68
Sensitivity	5.7 pg/mL
Expiration Date	Please enquire.

