

Product Datasheet

Porcine TGFb1 ELISA Kit (orb779572)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Porcine
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, platelet-poor plasma, tissue homogenates, cell culture supernates and other biological fluids
Uniprot ID	P07200
Sensitivity	11.9 pg/mL
Expiration Date	Please enquire.

