

Product Datasheet

Gallus GC ELISA Kit (orb779587)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Gallus
Range	23.44-1500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P68259
Sensitivity	7.59 pg/mL
Expiration Date	Please enquire.

