

Product Datasheet

Guinea pig IL10 ELISA Kit (orb779588)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Guinea pig
Range	3.13-200 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	Q9Z1Y5
Sensitivity	1.22 pg/mL
Expiration Date	Please enquire.

