

Product Datasheet

Canine IGF1 ELISA Kit (orb779647)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Canine
Range	78.13-5000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P33712
Sensitivity	31 pg/mL
Expiration Date	Please enquire.

