

## Product Datasheet

### Animal IgA ELISA Kit (orb779670)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Bovine
<b>Range</b>	0.79-50 µg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, saliva and other biological fluids
<b>Sensitivity</b>	0.34 ug/mL
<b>Expiration Date</b>	Please enquire.

