

Product Datasheet

Animal IGFBP3 ELISA Kit (orb779702)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Bovine |
| Range | 0.63-40 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids |
| Uniprot ID | P20959 |
| Sensitivity | 0.256 ng/mL |
| Expiration Date | Please enquire. |

