

Product Datasheet

Animal GAL8 ELISA Kit (orb779793)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Bovine
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	Q2YDD7
Sensitivity	0.118 ng/mL
Expiration Date	Please enquire.

