

Product Datasheet

Porcine FGF1 ELISA Kit (orb779819)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Porcine
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P20002
Sensitivity	6.1 pg/mL
Expiration Date	Please enquire.

