

Product Datasheet

Mouse NPFf ELISA Kit (orb779933)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	7.82-500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9WVA8
Sensitivity	13.5 pg/mL
Expiration Date	Please enquire.

