

Product Datasheet

Animal sIgA ELISA Kit (orb779944)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Bovine |
| Range | 0.79-50 µg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma, saliva and other biological fluids |
| Sensitivity | 0.34 ug/mL |
| Expiration Date | Please enquire. |

