

Product Datasheet

Human USP6NL ELISA Kit (orb780209)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 0.16-10 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Tissue homogenates and other biological fluids. |
| Uniprot ID | Q92738 |
| Sensitivity | 0.058 ng/mL |
| Expiration Date | Please enquire. |

