

Product Datasheet

Rat APOC2 ELISA Kit (orb780376)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Rat |
| Range | 0.32-20 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | G3V8D4 |
| Sensitivity | 0.137 ng/mL |
| Expiration Date | Please enquire. |

