

Product Datasheet

Human PYY3 ELISA Kit (orb780529)

| | |
|----------------------------|---|
| Description | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| Reactivity | Human |
| Range | 15.63-1000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma and other biological fluids |
| Uniprot ID | Q5JQD4 |
| Sensitivity | 4.82 pg/mL |
| Expiration Date | Please enquire. |

