

## Product Datasheet

### Guinea pig GT ELISA Kit (orb780972)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Guinea pig
<b>Range</b>	15.63-1000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<b>P06885</b>
<b>Sensitivity</b>	4.35 pg/mL
<b>Expiration Date</b>	Please enquire.

