

Product Datasheet

Animal GHRL ELISA Kit (orb781176)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Bovine
Range	156.25-10000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9BDJ6
Sensitivity	49.2 pg/mL
Expiration Date	Please enquire.

