

Product Datasheet

Rat POMC ELISA Kit (orb781458)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Rat
Range	156.25-10000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P01194
Sensitivity	52.5 pg/mL
Expiration Date	Please enquire.

