

Product Datasheet

Mouse UCHL1 ELISA Kit (orb781543)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.79-50 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9R0P9
Sensitivity	0.32 ng/mL
Expiration Date	Please enquire.

