

## Product Datasheet

### Guinea pig MCP1 ELISA Kit (orb781549)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Guinea pig
<b>Range</b>	15.63-1000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<a href="#">Q08782</a>
<b>Sensitivity</b>	6.8 pg/mL
<b>Expiration Date</b>	Please enquire.

